


# Organic solvents



System	: GC/FID	Solutes : 1. Methanol	26. 1,4-Dioxane
Column	: InertCap 17 0.25mmI.D. × 60m df=0.25µm	2. Ethyl ether	27. Toluene
Column Temp.	: 40°C (5min) → 4°C/min → 230°C (5min)	3. i-Propanol	28. Tetrachloroethylene
Carrier Gas	: He 130kPa	4. n-Hexane	29. Methyl-n-butyl ketone
Injection	: Split 100mL/min 1µL 250°C	5. Acetone	30. n-Butyl acetate
Detection	: FID Range 10 <sup>1</sup> 250°C	6. Carbon disulfide	31. p-Xylene
Sample	: Organic solvents Mixed evenly	7. Methyl acetate	32. m-Xylene
		8. Dichloromethane	33. Chlorobenzene
		9. trans-1,2-Dichloroethylene	34. N,N-Dimethyl formamide
		10. 2-Butanol	35. o-Xylene
		11. cis-1,2-Dichloroethylene	36. 1-Methylcyclohexanol
		12. Methyl ethyl keton	37. Cyclohexanol
		13. i-Butanol	38. Butyl cellosolve
		14. Ethyl acetate	39. Styrene
		15. Chloroform	40. Cellosolve acetate
		16. 1,1,1-Trichloroethane	41. Cyclohexanone
		17. Carbon Tetrachloride	42. 1,1,2,2-Tetrachloroethane
		18. Tetrahydrofuran	43. o-Dichlorobenzene
		19. Methylcellosolve	44. o-Cresol
		20. n-Butanol	45. m-Cresol
		21. 1,2-Dichloroethane	46. p-Cresol
		22. Trichloroethylene	
		23. n-Propyl acetate	
		24. i-Amyl alcohol	
		25. Ethyl cellosolve	