


Organic solvents



System : GC/FID
 Column : InertCap Pure-WAX
 0.25mmI.D. × 30m df=0.25µm
 Column Temp. : 40°C(5min)→5°C/min→200°C
 Carrier Gas : He 100kPa
 Injection : Split 1 : 80 0.1µl
 250°C
 Detection : FID Range 10⁰
 260°C
 Sample : Organic solvents

Solutes : 1. n-Hexane
 2. Cyclopentane
 3. Cyclohexane
 4. Cycloheptane
 5. Benzene
 6. Toluene
 7. Ethylbenzene
 8. o-Xylene
 9. α - Methylstyrene
 10. Cyclohexanol