


Analysis of Residual Solvent Components in Pharmaceutical Products – Using InertCap Pure-WAX

This application is an example of the analysis of components subject to the Guideline for Residual Solvents in Drugs.

Depending on the level of toxicity, this guideline classifies solvents as Class 1, Class 2, Class 3, and those for which no appropriate toxicity data are available. For these solvents, simultaneous analysis of the components was performed using GC/FID.

Example: measurement of mixed sample



Sample name and retention index

InertCap Pure WAX					
Component	RI	Component	RI	Component	RI
1. <i>n</i> -Pentane	500	25. 1,2-Dimethoxyethane	931	48. <i>p</i> -Xylene	1143
2. <i>n</i> -Hexane	600	26. 2-Propanol(Isopropyl alcohol)	931	49. 1-Butanol	1151
3. Diethyl ether	632	27. Dichloromethane	937	50. <i>m</i> -Xylene	1151
4. 2,2,4-Trimethylpentane	670	28. Ethanol	937	51. Nitromethane	1160
5. Diisopropyl ether	670	29. Methyl isopropyl ketone	937	52. 2-Methoxyethanol (Methyl cellosolve)	1186
6. <i>tert</i> -Butyl methyl ether	691	30. Benzene	948	53. Cumene	1186
7. <i>n</i> -Heptane	700	31. Propionaldehyde diethyl aceta	963	54. Pyridine	1191
8. 1,1-Dimethoxymethane	729	32. <i>n</i> -Propyl acetate	984	55. <i>o</i> -Xylene	1195
9. Cyclohexane	729	33. <i>cis</i> -1,2-Dichloroethylene	1001	56. 3-Methyl-1-butanol (Isoamyl alcohol)	1217
10. 1,1-Dichloroethylene	729	34. Trichloroethylene	1004	57. Chlorobenzene	1225
11. Methylcyclohexane	777	35. Acetonitrile	1012	58. 2-Ethoxyethanol	1231
12. Acetone	819	36. 4-Methyl-2-pentanone(MIBK)	1016	59. 1-Pentanol(Amyl alcohol)	1259
13. 2,2-Dimethoxypropane	830	37. Isobutyl acetate	1022	60. <i>N,N</i> -Dimethylformamide	1335
14. Methyl acetate	830	38. Chloroform	1031	61. Anisole	1354
15. Ethyl formate	830	39. 2-Butanol	1031	62. <i>N,N</i> -Dimethylacetamide	1413
16. Tetrahydrofuran	866	40. 1-Propanol	1044	63. Acetic acid	1448
17. <i>trans</i> -1,2-Dichloroethylen	866	41. Toluene	1049	64. 1,2,3,4-Tetrahydronaphthaler	1541
18. 2-Methyltetrahydrofuran	876	42. 1,4-Dioxane	1069	65. Dimethyl sulfoxide(DMSO)	1574
19. 1,1,1-Trichloroethane	891	43. 1,2-Dichloroethane	1074	66. Ethylene glycol	1631
20. Carbon tetrachloride	891	44. <i>n</i> -Butyl acetate	1083	67. <i>N</i> -methyl-2-pyrrolidone	1681
21. Ethyl acetate	897	45. 2-Hexanone(MBK)	1091	68. Formamide	1787
22. Methanol	903	46. 2-Methyl-1-propanol (Isobutyl alcohol)	1100		
23. 2-Butanone(MEK)	908				

Xylene used is a mixture of *m*-Xylene, *p*-Xylene, *o*-Xylene, and Ethylbenzene.

Retention indices are...

This value is based on the number of carbons in the straight-chain alkanes and is calculated using the retention time of each component.

In this application note, a temperature-rise analysis was made. The formula is shown below.

$$\text{Retention index} = 100 \times \frac{\text{TR} - \text{tR}(Z)}{\text{TR}(Z+1) - \text{tR}(Z)} + 100 \times Z$$

TR = retention time of the target compound

tR(Z) = retention time of straight-chain alkanes that precede the components of interest

TR(Z+1) = retention time of straight-chain alkanes emerging after the component of interest.

Z = number of carbons in the straight-chain alkane with hydrocarbon retention time tR(Z)

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-Chome
Shinjuku-ku, Tokyo,
163-1130, Japan
Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences B.V.

De Sleutel 9
5652 AS Eindhoven
The Netherlands
Phone: +31 (0)40 254 95 31
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences, Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
Phone: 310-265-4424
Fax: 310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

GL Sciences (ShangHai) Ltd.

Tower B, Room 2003,
Far East International Plaza,
NO,317 Xianxia Road,
Changning District.
Shanghai, China P.C. 200032
Phone: +86 (0)21-6278-2272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn

International Distributors

Visit our Website at www.glsciences.com/distributors

