


Analysis and Retention Index of 61 Components of Organic Solvents using Nitrogen Carrier Gas - Using InertCap 1

The retention index is a relatively representative index of the retention ratio of straight-chain alkanes and is used to study constituents based on the number of carbons in the molecule. It is one of the most useful pieces of information for qualitative analysis.

The retention index can be determined because in isothermal analysis the logarithm of the retention ratio for straight-chain alkanes is linearly related to the number of carbons, and the retention ratio is also linear to the number of carbons in thermal rise analysis.

In this report, InertCap 1 was used to determine the isothermal and temperature-rise retention indices of 61 organic solvents using nitrogen as the carrier gas.

Example: Measurement of standards



Conditions

System : GC - FID
Column : InertCap 1
 0.25 mm I.D. x 60 m df = 0.25 μ m
Col. Temp. : 40 $^{\circ}$ C - 5 $^{\circ}$ C/min - 220 $^{\circ}$ C
Carrier Gas : N₂ 90 kPa
Injection : Split 1:50
 240 $^{\circ}$ C
Detection : FID Range 10⁰
 240 $^{\circ}$ C
Sample Size : Mixed evenly
 0.2 μ L

Chromatographic conditions described above.

For isothermal analysis, adjust the pressure so that the linear velocity is constant.

Retention index in the temperature-rise analysis

| Peak No. | Component | Retention index | Retention time | Peak No. | Component | Retention index | Retention time |
|----------|--|-----------------|----------------|----------|--|-----------------|----------------|
| 1 | Methanol | 397 | 7.674 | 32 | 3-Methyl-1-butanol (Isoamyl alcohol) | 717 | 12.567 |
| 2 | Ethanol | 405 | 7.965 | 33 | 4-Methyl-2-pentanone (MIBK) | 719 | 12.634 |
| 3 | Acetonitrile | 428 | 8.084 | 34 | <i>N,N</i> -Dimethylformamide | 739 | 13.268 |
| 4 | Acetone | 448 | 8.183 | 35 | 1-Pentanol (Amyl alcohol) | 747 | 13.546 |
| 5 | 2-Propanol (Isopropyl alcohol) | 462 | 8.257 | 36 | Isobutyl acetate | 754 | 13.761 |
| 6 | Ethyl ether | 492 | 8.409 | 37 | Toluene | 755 | 13.782 |
| 7 | Carbon disulfide | 509 | 8.568 | 38 | 2-Hexanone (MBK) | 766 | 14.134 |
| 8 | <i>Tert</i> -Butanol | 506 | 8.528 | 39 | <i>N</i> -Butyl acetate | 795 | 15.071 |
| 9 | Methyl acetate | 508 | 8.554 | 40 | Tetrachloroethylene | 803 | 15.359 |
| 10 | Dichloromethane | 511 | 8.596 | 41 | <i>N,N</i> -Dimethylacetamide | 829 | 16.308 |
| 11 | 1-Propanol | 533 | 8.878 | 42 | Chlorobenzene | 832 | 16.441 |
| 12 | <i>Trans</i> -1,2-Dichloroethylene | 546 | 9.045 | 43 | Ethylbenzene | 851 | 17.128 |
| 13 | Methyl ethyl ketone | 568 | 9.327 | 44 | Isopentyl acetate (Isoamyl acetate) | 858 | 17.381 |
| 14 | 2-Butanol | 580 | 9.491 | 45 | <i>M</i> -Xylene | 859 | 17.448 |
| 15 | <i>cis</i> -1,2-Dichloroethylene | 589 | 9.602 | 46 | <i>P</i> -Xylene | 860 | 17.488 |
| 16 | Ethyl acetate | 596 | 9.697 | 47 | Cyclohexanone | 866 | 17.682 |
| 17 | <i>n</i> -Hexane | 600 | 9.746 | 48 | Cyclohexanol | 866 | 17.682 |
| 18 | Chloroform | 602 | 9.786 | 49 | 2-Ethoxyethyl acetate (Cellosolve acetate) | 877 | 18.12 |
| 19 | 2-Methyl-1-propanol (Isobutyl alcohol) | 610 | 9.975 | 50 | Styrene | 878 | 18.125 |
| 20 | Tetrahydrofuran | 613 | 10.037 | 51 | 1-Methylcyclohexanol | 879 | 18.177 |
| 21 | 2-Methoxyethanol (Methyl cellosolve) | 614 | 10.060 | 52 | 1,1,2,2-Tetrachloroethane | 881 | 18.269 |
| 22 | 1,2-Dichloroethane | 623 | 10.281 | 53 | <i>O</i> -Xylene | 883 | 18.332 |
| 23 | 1,1,1-Trichloroethane | 632 | 10.478 | 54 | 2-Butoxyethanol (Butyl cellosolve) | 888 | 18.528 |
| 24 | Isopropyl acetate | 641 | 10.683 | 55 | <i>N</i> -Pentyl acetate | 894 | 18.742 |
| 25 | 1-Butanol | 642 | 10.705 | 56 | 4-Methylcyclohexanone | 931 | 20.171 |
| 26 | Benzene | 654 | 10.975 | 57 | Phenol | 953 | 21.037 |
| 27 | Carbon tetrachloride | 653 | 10.948 | 58 | 1,2-Dichlorobenzene | 1018 | 23.582 |
| 28 | 1,4-Dioxane | 685 | 11.684 | 59 | <i>O</i> -Cresol | 1025 | 23.868 |
| 29 | Trichloroethylene | 687 | 11.723 | 60 | <i>P</i> -Cresol | 1048 | 24.752 |
| 30 | 2-Ethoxyethanol (Cellosolve) | 693 | 11.877 | 61 | <i>M</i> -Cresol | 1049 | 24.791 |
| 31 | <i>n</i> -Propyl acetate | 695 | 11.915 | | | | |

* Retention time in minutes

* Components with a symmetry factor of 1.5 or higher are highlighted in red.

In the case of temperature programming...

Because the retention ratio of straight-chain alkanes is linearly related to the number of carbons, the retention index is given by the following equation.

$$\text{Retention index } I = 100 \times \frac{\text{TR} - \text{tR}(Z)}{\text{TR}(Z+1) - \text{tR}(Z)} + 100 \times Z$$


TR = retention time of the target component

tR(Z) = retention time of straight-chain alkanes that precede the components of interest

tR(Z+1) = retention time of straight-chain alkanes emerging after the components of interest

Z = number of carbons in straight-chain alkanes with a retention time tR(Z)

Linear alkane (Z)



Retention index in isothermal analysis-1

| Peak No. (gradient temp.) | Component | 40 °C | | 80 °C | | 120 °C | | 160 °C | |
|------------------------------|---|-----------------|----------------|-----------------|----------------|-----------------|----------------|-----------------|----------------|
| | | Retention index | Retention time | Retention index | Retention time | Retention index | Retention time | Retention index | Retention time |
| 1 | Methanol | 378 | 7.619 | 354 | 7.195 | 317 | 7.745 | 387 | 7.502 |
| 2 | Ethanol | 448 | 8.037 | 418 | 7.318 | 415 | 7.815 | 447 | 7.539 |
| 3 | Acetonitrile | 460 | 8.141 | 445 | 7.383 | 447 | 7.850 | 476 | 7.563 |
| 4 | Acetone | 480 | 8.333 | 461 | 7.428 | 457 | 7.862 | 485 | 7.572 |
| 5 | 2-Propanol (Isopropyl alcohol) | 491 | 8.462 | 468 | 7.449 | 461 | 7.867 | 484 | 7.571 |
| 6 | Ethyl ether | 509 | 8.695 | 491 | 7.529 | 489 | 7.909 | 498 | 7.585 |
| 7 | Carbon disulfide | 509 | 8.687 | 491 | 7.529 | 491 | 7.912 | 493 | 7.580 |
| 8 | <i>Tert-Butanol</i> | 520 | 8.869 | 503 | 7.577 | 506 | 7.937 | 515 | 7.600 |
| 9 | Methyl acetate | 524 | 8.930 | 507 | 7.592 | 508 | 7.941 | 513 | 7.598 |
| 10 | Dichloromethane | 526 | 8.967 | 515 | 7.627 | 519 | 7.962 | 534 | 7.617 |
| 11 | 1-Propanol | 552 | 9.496 | 534 | 7.720 | 527 | 7.977 | 541 | 7.624 |
| 12 | <i>Trans-1,2-Dichloroethylene</i> | 562 | 9.725 | 555 | 7.838 | 556 | 8.043 | 575 | 7.661 |
| 13 | Methyl ethyl ketone | 581 | 10.249 | 572 | 7.953 | 571 | 8.081 | 588 | 7.677 |
| 14 | 2-Butanol | 592 | 10.585 | 580 | 8.008 | 582 | 8.112 | 590 | 7.680 |
| 15 | <i>Cis-1,2-Dichloroethylene</i> | 597 | 10.756 | 593 | 8.105 | 594 | 8.149 | 613 | 7.710 |
| 16 | Ethyl acetate | 604 | 11.021 | 592 | 8.097 | 587 | 8.128 | 592 | 7.683 |
| 17 | <i>N-Hexane</i> | 606 | 11.067 | 598 | 8.152 | 598 | 8.162 | 608 | 7.703 |
| 18 | Chloroform | 606 | 11.092 | 602 | 8.188 | 610 | 8.200 | 622 | 7.723 |
| 19 | 2-Methyl-1-propanol (Isobutyl alcohol) | 620 | 11.633 | 606 | 8.220 | 602 | 8.175 | 613 | 7.711 |
| 20 | Tetrahydrofuran | 619 | 11.589 | 617 | 8.323 | 624 | 8.249 | 640 | 7.751 |
| 21 | 2-Methoxyethanol (Methyl cellosolve) | 620 | 11.633 | 613 | 8.283 | 615 | 8.218 | 629 | 7.733 |
| 22 | 1,2-Dichloroethane | 630 | 12.099 | 630 | 8.455 | 638 | 8.301 | 654 | 7.773 |
| 23 | 1,1,1-Trichloroethane | 634 | 12.277 | 640 | 8.567 | 650 | 8.348 | 666 | 7.794 |
| 24 | Isopropyl acetate | 649 | 13.063 | 642 | 8.584 | 636 | 8.291 | 644 | 7.756 |
| 25 | 1-Butanol | 649 | 13.069 | 642 | 8.580 | 638 | 8.302 | 645 | 7.759 |
| 26 | Benzene | 649 | 13.061 | 657 | 8.766 | 666 | 8.423 | 684 | 7.828 |
| 27 | Carbon tetrachloride | 655 | 13.371 | 662 | 8.838 | 673 | 8.452 | 690 | 7.840 |
| 28 | 1,4-Dioxane | 685 | 15.442 | 689 | 9.233 | 695 | 8.567 | 711 | 7.883 |
| 29 | Trichloroethylene | 686 | 15.488 | 692 | 9.280 | 699 | 8.594 | 716 | 7.894 |
| 30 | 2-Ethoxyethanol (Cellosolve) | 694 | 16.187 | 693 | 9.298 | 691 | 8.548 | 704 | 7.868 |
| 31 | <i>N-Propyl acetate</i> | 699 | 16.563 | 693 | 9.302 | 689 | 8.536 | 696 | 7.851 |

* Retention time in minutes

* Components with a symmetry factor of 1.5 or higher are highlighted in red.

In the case of isothermal analysis...

Because the logarithm of the retention ratio of straight-chain alkanes is linearly related to the number of carbons, the retention index is given by the following equation.

$$\text{Retention index } I = 100 \times \frac{\log t'R - \log t'R(Z)}{\log t'R(Z+1) - \log t'R(Z)} + 100 \times Z$$

$t'R$ = retention time of the target component

$t'R(Z)$ = retention time of straight-chain alkanes that precede the components of interest

$t'R(Z+1)$ = retention time of straight-chain alkanes emerging after the components of interest

Z = number of carbons in straight-chain alkanes with a retention time $t'R(Z)$

$t'R$ = corrected retention time $t'R = t'R - t'O$

$t'O$ = hold-up time (elution time of non-retentive components)

Retention index in isothermal analysis-2

| Peak No. (gradient temp.) | Component | 40 °C | | 80 °C | | 120 °C | | 160 °C | |
|------------------------------|---|-----------------|----------------|-----------------|----------------|-----------------|----------------|-----------------|----------------|
| | | Retention index | Retention time | Retention index | Retention time | Retention index | Retention time | Retention index | Retention time |
| 32 | 3-Methyl-1-butanol (Isoamyl alcohol) | 722 | 18.851 | 715 | 9.712 | 713 | 8.673 | 722 | 7.908 |
| 33 | 4-Methyl-2-pentanone (MIBK) | 721 | 18.697 | 721 | 9.827 | 725 | 8.748 | 736 | 7.941 |
| 34 | <i>N,N</i> -Dimethylformamide | 739 | 20.861 | 742 | 10.301 | 748 | 8.905 | 765 | 8.019 |
| 35 | 1-Pentanol(Amyl alcohol) | 754 | 22.924 | 747 | 10.423 | 745 | 8.884 | 752 | 7.983 |
| 36 | Isobutyl acetate | 760 | 23.753 | 755 | 10.623 | 753 | 8.946 | 758 | 8.000 |
| 37 | Toluene | 752 | 22.648 | 760 | 10.790 | 771 | 9.091 | 787 | 8.084 |
| 38 | 2-Hexanone(MBK) | 768 | 25.064 | 768 | 10.988 | 769 | 9.077 | 781 | 8.067 |
| 39 | <i>n</i> -Butyl acetate | 799 | 31.041 | 794 | 11.813 | 789 | 9.259 | 795 | 8.111 |
| 40 | Tetrachloroethylene | 795 | 30.203 | 806 | 12.255 | 819 | 9.562 | 836 | 8.259 |
| 41 | <i>N,N</i> -Dimethylacetamide | 827 | 38.136 | 830 | 13.123 | 838 | 9.787 | 851 | 8.318 |
| 42 | Chlorobenzene | 823 | 36.940 | 834 | 13.332 | 849 | 9.931 | 867 | 8.388 |
| 43 | Ethylbenzene | 843 | 43.012 | 852 | 14.222 | 864 | 10.140 | 878 | 8.442 |
| 44 | Isopentyl acetate (Isoamyl acetate) | 862 | 49.713 | 858 | 14.520 | 856 | 10.028 | 859 | 8.352 |
| 45 | <i>m</i> -Xylene | 852 | 45.945 | 861 | 14.677 | 870 | 10.234 | 882 | 8.463 |
| 46 | <i>p</i> -Xylene | 853 | 46.363 | 862 | 14.732 | 871 | 10.242 | 885 | 8.478 |
| 47 | Cyclohexanone | 854 | 46.568 | 866 | 15.008 | 882 | 10.421 | 902 | 8.562 |
| 48 | Cyclohexanol | 860 | 49.053 | 866 | 14.958 | 875 | 10.309 | 890 | 8.501 |
| 49 | 2-Ethoxyethyl acetate (Cellosolve acetate) | 886 | 60.426 | 878 | 15.695 | 871 | 10.253 | 871 | 8.408 |
| 50 | Styrene | 868 | 52.032 | 878 | 15.709 | 889 | 10.542 | 905 | 8.578 |
| 51 | 1-Methylcyclohexanol | 870 | 52.877 | 879 | 15.772 | 890 | 10.559 | 908 | 8.597 |
| 52 | 1,1,2,2-Tetrachloroethane | 873 | 54.268 | 881 | 15.920 | 892 | 10.589 | 910 | 8.607 |
| 53 | <i>o</i> -Xylene | 872 | 53.872 | 883 | 16.058 | 895 | 10.648 | 912 | 8.617 |
| 54 | 2-Butoxyethanol (Butyl cellosolve) | 887 | 61.049 | 888 | 16.368 | 890 | 10.554 | 897 | 8.537 |
| 55 | <i>n</i> -Pentyl acetate | 899 | 66.827 | 894 | 16.809 | 891 | 10.570 | 891 | 8.504 |
| 56 | 4-Methylcyclohexanone | — | — | 928 | 19.542 | 945 | 11.648 | 967 | 8.966 |
| 57 | Phenol | — | — | 955 | 22.227 | 946 | 11.681 | 948 | 8.841 |
| 58 | 1,2-Dichlorobenzene | — | — | 1013 | 30.331 | 1031 | 14.255 | 1052 | 9.724 |
| 59 | <i>o</i> -Cresol | — | — | 1028 | 33.057 | 1024 | 14.001 | 1031 | 9.509 |
| 60 | <i>p</i> -Cresol | — | — | 1049 | 37.332 | 1044 | 14.762 | 1047 | 9.665 |
| 61 | <i>m</i> -Cresol | — | — | 1050 | 37.593 | 1045 | 14.808 | 1047 | 9.672 |

* Retention time in minutes

* Components with a symmetry factor of 1.5 or higher are highlighted in red.

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-Chome
Shinjuku-ku, Tokyo,
163-1130, Japan
Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences B.V.

De Sleutel 9
5652 AS Eindhoven
The Netherlands
Phone: +31 (0)40 254 95 31
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences, Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
Phone: 310-265-4424
Fax: 310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

GL Sciences (ShangHai) Ltd.

Tower B, Room 2003,
Far East International Plaza,
NO,317 Xianxia Road,
Changning District.
Shanghai, China P.C. 200032
Phone: +86 (0)21-6278-2272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn

International Distributors

Visit our Website at www.glsciences.com/distributors

