

Using MonoTrap RGC18 TD and HandyTD TD265, volatile components in Japanese sake were screened and analyzed, with differences in the analytical results compared between the two different collection methods.

MonoTrap RGC18 TD is used to collect volatile constituents. HandyTD TD265 is used to introduce the volatile compounds collected on the trapping agent into a GC by thermal desorption. Volatile components were collected by a dip-shaking method (a method in which a MonoTrap is added to sake, with the volatile components collected by shaking with a constant-temperature shaker at 60 °C). The effect of salting-out was confirmed by comparing a method of sample collection without salting out, with a method with sodium chloride added to the sample.

Salting-out improves the sensitivity for all components with detection of many flavor compounds in Japanese sake. In particular, Isovaleraldehyde diethyl acetal, a fragrance substance produced by the deterioration of sake, was detected by salting-out.

Pretreatment procedure

Japanese sake (rice wine)

Place 15 mL of Japanese Sake in 40 mL vials

GC/MS Conditions

System

: GC - MS - Thermal Desorption (HandyTD TD265)

Column

: InertCap Pure-WAX
0.25 mm I.D. x 60 m, df = 0.5 μm

Col. Cat. No.

: 1010-68164

Col.Temp.

: 40 °C (5 min) - 10 °C/min - 250 °C

Carrier Gas

: He, 1 mL/min (constant flow)

GC Inlet

: 250 °C Split 10:1

Detection

: MS Scan (m/z 30-350)

Collection (HS)

MonoTrap RGC18 TD × 1

60 °C 1 hour collection



Sample 1
With salting-out



Sample 2
No salting out

HandyTD/GC/MS



HandyTD Conditions


Desorb Temp.

: Room temperature - 45 °C/sec -


Pre Desorb Press.

: 200 °C (5 min)

: 140 kPa



Analysis results of sake with salt analysis



* Standard samples are not used for qualitative analysis.
Results from a library search.

- | | | |
|------------------------------------|----------------------------------|------------------------------|
| 1. Ethyl acetate | 17. Propylene glycol ethyl ether | 33. Benzeneacetaldehyde |
| 2. Methanol | 18. Acetic acid | 34. Diethyl succinate |
| 3. Ethanol | 19. Furfural | 35. Dodecanal |
| 4. Isobutyl acetate | 20. Formic acid | 36. Phenylethyl formate |
| 5. 1-Propanol | 21. Ethyl hydroxybutyrate | 37. Caproic acid |
| 6. Ethyl butyrate | 22. 2,3-Butanediol | 38. Phenethyl acetate |
| 7. Isovaleraldehyde diethyl acetal | 23. Ethyl 2-hydroxycaproate | 39. Phenylethyl alcohol |
| 8. Isobutanol | 24. 1-Octanol | 40. Tetradecanal |
| 9. Isoamyl acetate | 25. Benzaldehyde | 41. 1-Dodecanol |
| 10. Butanol | 26. 2,3-Butanediol | 42. Octanoic acid |
| 11. 2-Methyl-1-butanol | 27. 2-Oxopentanedioic acid | 43. Dihydroxyacetone |
| 12. Isoamyl alcohol | 28. 5-Methyl-furfural | 44. Diisobutyl adipate |
| 13. Ethyl caproate | 29. 2-Cyclopentene-1,4-dione | 45. Myristyl alcohol |
| 14. Acetoin | 30. Butanoic acid | 46. Hydroxydihydromaltol |
| 15. Acetol | 31. Furfuryl alcohol | 47. Glycerin |
| 16. Ethyl lactate | 32. Butyrolactone | 48. 5-Hydroxymethyl furfural |

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-Chome
Shinjuku-ku, Tokyo,
163-1130, Japan
Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences B.V.

De Sleutel 9
5652 AS Eindhoven
The Netherlands
Phone: +31 (0)40 254 95 31
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences, Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
Phone: 310-265-4424
Fax: 310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

GL Sciences (ShangHai) Ltd.

Tower B, Room 2003,
Far East International Plaza,
NO.317 Xianxia Road,
Changning District.
Shanghai, China P.C. 200032
Phone: +86 (0)21-6278-2272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn

International Distributors

Visit our Website at www.glsciences.com/distributors