

Enrichment Analysis of Volatile Components from Hairbrush — Sampling Comparisons Between MonoTrap and Sampling Tube


This technical note shows the differences in analysis by using two different sampling medias: MonoTrap RGC18 TD and Sampling tube (TenaxTA 150 g).

The method used this time showed better sensitivity for low boiling point components, such as Toluene and Xylene, when using sampling tube. On the other hand, MonoTrap showed better sensitivity for hydrocarbons, such as Hexadecane and the medium to high-boiling components (Methylbenzothiazole and Phenanthrene). Especially when placed on the hairbrush, the MonoTrap sensitivity was higher.

When using sampling tube, it is possible to sampling the entire amount of volatilized compounds in Tedlar Bag. Therefore, it could deliver a good sensitivity for easily volatilized compounds. On the contrary, as MonoTrap can be put very close to the sample, it can be used to achieve higher sensitivity of those compounds with low volatility.

Sampling Preparation

<Sampling with MonoTrap>




Place RGC18 (i)* on the hairbrush and put into Tedlar bag.
Place another RGC18 (ii)* at a distance from the hair brush.
Remove the air from the bag with an air pump.
Add nitrogen 1 L.

Collect headspace gas overnight at room temperature




* RGC18 (i) is in contact with the hairbrushes
RGC18 (ii) is not in contact with hairbrushes

<Sampling with sampling tube>



Place a hairbrush in a tedra bag
Remove air from the pump
Add 1 L

Allow to stand at room temperature overnight to volatilize the volatile components
Connect the collection tube to the sampling pump
Connect Tedra Bags behind Collection Tubes
Collect 1 L (100 mL/min) with sampling pumps



※ Volatile constituents in Tedra bags sampled
Be aspirated by pumps and collected in the TenaxTA

System Requirements


GC/MS Conditions

System	: Thermal Desorption-GC/MS (HandyTD TD265)
Column	: InertCap Pure-WAX (0.25 mm I.D. × 60 m, df = 0.5 μm)
Col. Cat. No.	: 1010-68164
Col. Temp.	: 40 °C(5 min) - 5 °C/min - 250 °C
Carrier Gas	: He, 1 mL/min (constant flow)
GC Inlet	: 250 °C Split 10:1
Detection	: MS Scan (<i>m/z</i> 30-350)

HandyTD Conditions


Desorb Temp.	: Room temperature(0.1 min) - 45 °C/sec- 250 °C(5 min) ⇒ MonoTrap
Pre Desorb Press.	: Room temperature(1 min) - 45 °C/sec- 270 °C(5 min) ⇒ Tenax TA : 140 kPa

Comparison of Operating Blanks




* 1, Dimethylacetamide; 2, Phenol, ingredients from tetra bags.

Comparison of Analysis Result



Analytical results (enlarged) of RGC18 (a) and TenaxTA.



No qualitative tests have been performed on standard samples.
Library search results.

- | | |
|-----------------------------|---|
| 1. Pinene | 16. Heptadecane |
| 2. Toluene | 17. Teramethyl Hexadecane |
| 3. Pinene | 18. Octadecane |
| 4. Ethylbenzene | 19. Benzyl Alcohol |
| 5. Xylene | 20. Propanoic acid, 2-methyl-, 1-(1,1-dimethylethyl)-2-methyl-1,3-propanediyl ester |
| 6. Xylene | 21. Butylated Hydroxytoluene |
| 7. Xylene | 22. Methylbenzothiazole |
| 8. Cyclohexanone | 23. Benzothiazole |
| 9. Dimethylformamide | 24. Phenol* |
| 10. Methylstyrene | 25. Acetone anil |
| 11. Acetyl dimethylcarbinol | 26. Caprolactam |
| 12. Butoxyethanol | 27. Fluorene |
| 13. Dimethylacetamide* | 28. Phenanthrene |
| 14. Pentadecane | |
| 15. Hexadecane | |

* This is an ingredient derived from a tetra bag.

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan

22-1 Nishishinjuku 6-chome
Shinjuku-ku, Tokyo
163-1130, Japan

Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
USA

Phone: +1-310-265-4424
Fax: +1-310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

GL Sciences B.V.

Dillenburgstraat 7C
5652AM, Eindhoven
The Netherlands

Phone: +31-40-254-9531
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (Shanghai) Limited

Tower B, Room 2003
Far East International Plaza
No.317 Xianxia Road, Changning District
Shanghai, China 200051

Phone: +86-21-62782272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn



International Distributors

Visit our Website at www.glsciences.com/distributors