


Analysis of Sugars(detected by Post Column method)



System	: GL-7400 system	Solutes	: 1. Rhamnose (100 µg/mL)
Column	: Asahipak NH2P-50 4E (250 x 4.6 mm I.D.)		2. Ribose (100 µg/mL)
Column Cat. No.	: 5055-12611		3. Fucose (100 µg/mL)
Eluent	:		4. Xylose (100 µg/mL)
Flow rate	: 1.0 mL/min.		5. Arabinose (100 µg/mL)
Col. Temp.	: 40 °C		6. Sorbose (100 µg/mL)
Detection	: FL Ex 330 nm Em 470 nm (GL-7453 FL Detector)		7. Fructose (100 µg/mL)
Injection Vol.	: 10 µL		8. Mannose (100 µg/mL)
Sample	: Sugars		9. Glucose (100 µg/mL)
			10. Galactose (100 µg/mL)
			11. Sucrose (100 µg/mL)
			12. Maltose (100 µg/mL)
			13. Lactose (100 µg/mL)
			14. Trehalose (100 µg/mL)
			15. Kestose (100 µg/mL)
			16. Raffinose (100 µg/mL)
			17. Stachyose (100 µg/mL)