


Analysis of Organic acids(detected by Post-Column method)

Data No. LA539-0000



System	: GL-7400 system	Solutes	: 1. Phosphoric acid (1.0 mg/mL)
Column	: Inertsil CX (5 µm, 250 × 4.6 mm I.D.) × 2		2. Citric acid (1.0 mg/mL)
Column Cat. No.	: 5020-07146		3. Malic acid (1.0 mg/mL)
Eluent	: 0.1 % TFA		4. Succinic acid (1.0 mg/mL)
Flow rate	: 1.0 mL/min.		5. Lactic acid (1.0 mg/mL)
Col. Temp.	: 55 °C		6. Acetic acid (1.0 mg/mL)
Detection	: VIS 440 nm (GL-7450 UV Detector)		7. Pyroglutamic acid (1.0 mg/mL)
Injection Vol.	: 10 µL		
Sample	: Organic acids		