


## Analysis of Glycyrrhiza Standard (Under the Condition of Japanese Pharmacopoeia)

Data No. LA550-0698



System	:GL-7400 system	Solutes	: 1. Glycyrrhizic acid dipotassium salt (0.27 mg/mL)
Column	:Inertsil ODS-SP (5 $\mu$ m, 150 $\times$ 4.6 mm I.D.)		2. Propyl p-hydroxybenzoate (0.05 mg/mL)
Column Cat. No.:			
Eluent	: A) CH <sub>3</sub> CN B) 2.1 % CH <sub>3</sub> COOH A/B = 40/60 ,v/v		
Flow rate	:0.6 mL/min		
Col. Temp.	:30 °C		
Detection	:PDA 254 nm (GL-7452A PDA detector)		
Injection Vol.	:20 $\mu$ L		
Sample	:Glycyrrhiza Standard		