


## Analysis of Iminoctadine acetate by Post-column derivatization

Data No. LA570-0696



System :GL-7400 system Solutes :1. Iminoctadine acetate (12  $\mu$ g/L)

Column :Inertsil ODS-3 (5  $\mu$ m, 150 x 4.6 mm I.D.)

Column Cat. No.:

Eluent : A) CH<sub>3</sub>CN  
B) 0.1 M NaClO<sub>4</sub>, 10 mM NaOH, 24 mM Lactic acid  
A/B = 5/17,v/v

Flow rate :1.2 mL/min

Col. Temp. :40 °C

Detection :FL Ex 395 nm Em 500 nm (GL-7453A FL Detector)

Injection Vol. :200  $\mu$ L

Sample :