


# InertSearch™ for LC

Inertsil® Applications

## Analysis of Polycyclic aromatic hydrocarbons (PAHs) (Using UV Detector)

Data No. LA732-0871



### Conditions

**System** : GL-7400 HPLC system  
**Column** : Inertsil ODS-P HP (3  $\mu$  m, 100 x 3.0 mm I.D.)  
**Column Cat. No.** : 5020-87036  
**Eluent** : A) CH<sub>3</sub>OH  
B) H<sub>2</sub>O  
A/B = 80/20 - 2 min- 80/20 - 3 min - 100/0, v/v  
**Flow rate** : 0.6 mL/min  
**Col. Temp.** : 40 °C  
**Detection** : UV 254 nm (GL-7452A PDA Detector)  
**Injection Vol.** : 5  $\mu$  L  
**Sample** : PAHs

Reference [LA337](#), [LA607](#)

### Analyte:

1. Naphthalene	(5 mg/L)
2. Acenaphthylene	(5 mg/L)
3. 1-Methylnaphthalene	(5 mg/L)
4. Acenaphthene	(5 mg/L)
5. Fluorene	(5 mg/L)
6. Phenanthrene	(5 mg/L)
7. Anthracene	(5 mg/L)
8. Fluoranthene	(5 mg/L)
9. Pyrene	(5 mg/L)
10. Triphenylene	(5 mg/L)
11. Benzo-[a]-anthracene	(5 mg/L)
12. Chrysene	(5 mg/L)
13. Benzo-[b]-fluoranthene	(5 mg/L)
14. Benzo-[k]-fluoranthene	(5 mg/L)
15. Benzo-[a]-pyrene	(5 mg/L)
16. DiBenz-[a,h]-anthracene	(5 mg/L)
17. Benzo-[g,h,i]-perylene	(5 mg/L)
18. Indeno-[1,2,3-cd]pyrene	(5 mg/L)