


## Analysis of Clopidogrel sulfate

(Under the Condition of the Japanese Pharmacopoeia)



### Conditions

**System** : GL-7400 HPLC system

**Column** : Inertsil WP300 C18  
(5  $\mu$  m, 150 x 3.9 mm I.D.)

**Column Cat. No.** : 5020-

**Eluent** : A) 0.87 g/L IPCC-05 in H<sub>2</sub>O (pH 2.5, H<sub>3</sub>PO<sub>4</sub>)/CH<sub>3</sub>OH = 19/1, v/v  
(IPCC-05: Sodium 1-Pentanesulfonate)  
B) CH<sub>3</sub>CN/CH<sub>3</sub>OH = 19/1, v/v  
A/B = 60/40, v/v

**Flow Rate** : 0.8 mL/min

**Col. Temp.** : 30 °C

**Detection** : UV 220 nm (GL-7452 PDA Detector)

**Injection Vol.** : 10  $\mu$  L

**Sample** : Standard

### Analyte:

1. Clopidogrel sulfate

126 mg/L

Theoretical plates : 9,973 ( $\geq$  4,500)

Tailing factor : 0.91 ( $\leq$  2.0)

RSD of the  
peak area (%) (n=6) : 0.18 ( $\leq$  1.0)