


Analysis of Acesulfame potassium

(Under the Condition of the JSFA-VIII (Japanese Standards of Food Additives))



Conditions

System : GL7700 HPLC system
Column : InertSustain C18
 (5 μ m, 250 x 4.6 mm I.D.)
Column Cat. No. : 5020-07346
Eluent : A) CH₃CN
 B) 10mM Tetrabutylammonium hydrogensulfate
 A/B = 2/3, v/v
Flow Rate : 0.83 mL/min
Col. Temp. : 40 °C
Detection : UV 227 nm (UV7750 UV Detector)
Injection Vol. : 20 μ L
Sample : Standard

Analyte:

1. Acesulfame potassium 10 mg/L
 2. *p*-Hydroxybenzoic acid ethyl ester 20 mg/L

Resolution (1, 2) : 15.0 (\geq 4)