


Analysis of bromate (Detected by Post-Column method)



Conditions

System	: Chromaster	Analyte:	
Column	: InertSustain C18 (5 μ m, 250 x 4.6 mmI.D.)	1. Bromate anion	20 μ g/L as KBrO ₃
Column Cat. No.	: 5020-07346		
Eluent	: Solution A		
Flow Rate	: 0.9 mL/min		
Reaction Reagent	: Solution B		
Reaction Temp.	: 60 °C		
Reaction flow Rate	: 0.3 mL/min		
Col. Temp.	: 40 °C		
Detection	: VIS 450 nm (5430 Diode Array Detector)		
Injection Vol.	: 100 μ L		
Sample	: KBrO ₃ Standard solution		

Solution A : Dissolve 2.0 g of CH₃COOH and 45 g of 10% Tetra-*n*-butylammonium Hydroxide (TBAH) aq. in mixture of methanol : water (700 mL : 100 mL).

Adjust with TBAH to pH 5.0, and make up to 1000 mL by adding water.

Solution B : A) Dissolve 10.0 g of KBr in mixture of water : nitric acid (700 mL : 60mL)

B) Dissolve 500 mg *o*-dianisidine dihydrochloride in 200mL of methanol

Mix of A and B, and make up to 1000 mL by adding water.