

Analysis of Colorants by HPLC-PDA


Artificial colorants are widely used because of their low prices and useful chemical properties. For example, as compared with natural food dyes, many artificial colorants are more stable to heat, light, and oxygen. This feature is advantageous in quality control of the colorants. However, it has been pointed out that their excess or combinational use may cause toxic effect. In this note, a simultaneous determination method for 18 artificial colorants using PDA (photo-diode array) detector is described.

Conditions

Column : Inertsil ODS-3 (5 μ m, 150 x 4.6 mm I.D. Cat.No. 5020-01731)
Eluent : A) CH₃CN, B) 10 mM Na₂HPO₄ (pH 6.9) A/B = 10/90 – 50 min - 35/65
Flow rate : 1.0 mL/min
Col. Temp. : 40 °C
Detector : PDA
Injection Vol. : 10 μ L

PDA detector is highly recommended for simultaneous determination of various compounds. Even if multiple peaks are detected at same retention time by monitoring with a certain wavelength, it may be possible to detect each compound as a single peak by detecting at other wavelengths.

The chromatogram shown below is an example. Peak 17 and peak 18 were not separated at wavelength of 270 nm. However, with shifting the wavelength to 470 nm peak 17 was obtained as single peak, and peak 18 also can be determined by detection at 580 nm.

**Components**

Red numbers : dyes not permitted as food additives by the Food Sanitation Law of Japan

1. Tartrazine (FD&C Yellow No.5, Acid Yellow 23)	7.6mg/L	10. Ponceau SX (FD&C Red No.4, Scarlet GN)	5.3 mg/L
2. Amaranth (FD&C Red No.2, Acid Red 27)	3.8mg/L	11. Orange I (Ext. D&C Orange No.3, Acid Orange 20)	5.3 mg/L
3. Indigo Carmine (FD&C Blue No.2, Acid Blue W)	7.6mg/L	12. Fast Green FCF (FD&C Green No.3, Solid Green FCF)	
4. Ponceau 4R (Acid Red 18)	3.8mg/L	13. Brilliant Blue FCF (FD&C Blue No.1, Acid Blue 9)	3.0 mg/L
5. Sunset Yellow FCF (FD&C Yellow No. 6, Orange Yellow S)	5.3 mg/L	14. Ponceau 3R (FD&C Red No.1)	3.0mg/L
6. Naphthol Yellow S	7.6 mg/L	15. Erythrosine (FD&C Red No.3, Acid Red 51)	7.6mg/L
7. Uranine (D&C Yellow No.8, Fluorescein Sodium Salt)	3.8mg/L	16. Azure Blue VX (Acid Blue 1, Patent Blue V)	5.3mg/L
8. Allura Red AC (FD&C Red No.40)	5.3mg/L	17. Orange II (D&C Orange No.4, Acid Orange 7)	3.0mg/L
9. Ponceau R (D&C Red No.5, Acid Red 26)	7.6mg/L	18. Sulforhodamine B (Acid Red 52)	7.6mg/L
			3.0mg/L

Analysis of Powdered Soft Drink (powdered coke)

Pretreatment Procedure

Powdered coke 10 g


Dissolution

10 mL Water

Filtration

Membrane filter
(0.45 μ m)

HPLC



GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-Chome
Shinjuku-ku, Tokyo,
163-1130, Japan
Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences B.V.

De Sleutel 9
5652 AS Eindhoven
The Netherlands
Phone: +31 (0)40 254 95 31
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (ShangHai) Ltd.

Tower B, Room 2003,
Far East International Plaza,
NO,317 Xianxia Road,
Changning District.
Shanghai, China P.C. 200032
Phone: +86 (0)21-6278-2272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn

GL Sciences, Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
Phone: 310-265-4424
Fax: 310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

International Distributors

Visit our Website at:

<https://www.glsciences.com/company/distributor.html>