


The analysis of geniposide is presented in this application note. Geniposide is obtained from gardenia (hawthorn) and is used as a food additive derived from natural products. It is also considered to be a crude drug listed in the Japanese Pharmacopoeia.



Analysis was performed under the conditions described in the draft formulary for food additives, 8th edition. As a result, it was possible to easily detect geniposide without any influence from contamination.

(K.Suzuki)

Structural formula



Structures are created using Chemistry 4-D Draw which is Provided by ChemInnovation Software, Inc.

Example: Measurement of standard**1. Geniposide (1.0 mg/L)****HPLC conditions**

Column : Inertsil ODS-SP
(5 μ m, 150 x 3.0 mm I.D.)


Eluent : A) CH₃ CN
B) H₂O
A/B = 15/85, v/v
(gradient mixer)


Flow rate : 0.4 mL/min

Column temperature : 40 °C

Detected : PDA 240 nm

Injection volume : 20 μ L



Examples: Analysis of gardenia extracts

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-Chome
Shinjuku-ku, Tokyo,
163-1130, Japan
Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences B.V.

De Sleutel 9
5652 AS Eindhoven
The Netherlands
Phone: +31 (0)40 254 95 31
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (ShangHai) Ltd.

Tower B, Room 2003,
Far East International Plaza,
NO,317 Xianxia Road,
Changning District.
Shanghai, China P.C. 200032
Phone: +86 (0)21-6278-2272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn

GL Sciences, Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
Phone: 310-265-4424
Fax: 310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

International Distributors

Visit our Website at www.glsciences.com/distributors