


Analysis of Stevioside by HPLC with Photodiode Array Detection


Stevioside is a widely known natural sweetener extracted from the stevia of the family Asteraceae. The major glycoside is rebaudioside A.

The separated peaks were successfully detected using a PDA detector and an aminopropyl column.

(K.Suzuki)



Example: Measurement of standard



HPLC conditions

Column : Inertsil NH2
(5 μ m, 150 x 4.6 mm I.D.)


Eluent : A) CH₃ CN
B) H₂O
A/B = 80/20, v/v
(gradient mixer)

Flow rate : 1.0 mL/min

Column temperature : 40 °C


Detected : PDA 210 nm

Injection volume : 20 μ L



Examples: Analysis of soft drinks

The analysis was performed by filtering the samples of soft drinks through a 0.45 μm filter and diluting them equally with acetonitrile.



GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-Chome
Shinjuku-ku, Tokyo,
163-1130, Japan
Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences B.V.

De Sleutel 9
5652 AS Eindhoven
The Netherlands
Phone: +31 (0)40 254 95 31
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (ShangHai) Ltd.

Tower B, Room 2003,
Far East International Plaza,
NO,317 Xianxia Road,
Changning District.
Shanghai, China P.C. 200032
Phone: +86 (0)21-6278-2272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn

GL Sciences, Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
Phone: 310-265-4424
Fax: 310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

International Distributors

Visit our Website at www.glsciences.com/distributors