



Uric acid is the final metabolite of the purines in the body. However, when uric acid is present in excess in the body, gout, arthritis, and gouty kidney are known to develop. Serum levels of uric acid (or uric acid levels), which are used to diagnose abnormalities in purine metabolism and kidney dysfunction.

The uricase-UV method is a practical reference method for the measurement of uric acid. Although many methods have been proposed, an HPLC method with UV detection provides accurate and sensitive measurement.

This application presents an analytical method using Inertsil ODS-4 analytical column based on the recommended method (Clinical Chemistry 22:300-307, 1993). Analysis of serum samples revealed a well-shaped uric acid peak within 10 minutes.

(C. Aoyama, K. Suzuki)

Example: Measurement of standard



Calibration curve for uric acid


HPLC conditions

Guard column	: Cartridge guard column E Inertsil ODS-4 (5 μm, 10 x 4.0 mm I.D.)
Column	: Inertsil ODS-4 (5 μm, 250 x 4.6 mm I.D.)
Eluent	: A) 74 mM phosphate Buffer (pH 2.2)* B) CH ₃ OH A / B = 98 / 2, v/v
Flow rate	: 1.0 mL/min
Column temperature	: 40 °C
Detector	: UV 284 nm
Injection volume	: 50 μL


* 74 mM phosphate buffer (pH 2.2):
0.2 M disodium hydrogen phosphate (aqueous)
70 mL
and 0.5 M phosphoric acid (aqueous) 120 mL
was mixed with Ultrapure water and made up to a
total volume of 1 L.

Example of sample measurement

Example of pretreatment




Structural Formula



Uric acid

Structures are created using Chemistry 4-D Draw which is provided by ChemInnovation Software, Inc.

Measurement of uric acid in human control serum



GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-Chome
Shinjuku-ku, Tokyo,
163-1130, Japan
Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences B.V.

De Sleutel 9
5652 AS Eindhoven
The Netherlands
Phone: +31 (0)40 254 95 31
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (ShangHai) Ltd.

Tower B, Room 2003,
Far East International Plaza,
NO,317 Xianxia Road,
Changning District.
Shanghai, China P.C. 200032
Phone: +86 (0)21-6278-2272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn

GL Sciences, Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
Phone: 310-265-4424
Fax: 310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

International Distributors

Visit our Website at www.glsciences.com/distributors