



Clopidogrel is an antiplatelet agent. It is used for treatment of thrombosis, such as stroke and myocardial infarction. Aspirin, which is famous as a non-steroidal anti-inflammatory drug (NSAID), also has an antiplatelet effect.


Advantages and safety of combination drug of clopidogrel and aspirin have been studied and discussed. In this note, clopidogrel and aspirin were successfully separated by using Inertsil ODS-4 with superior peak shapes.

Chemical Structure



Structures are created using Chemistry 4-D Draw which is provided by ChemInnovation Software, Inc.

Standard Solution




Conditions

Column	: Inertsil ODS-4 HP (3 μm, 50 x 4.6 mm I.D.)
Column Cat. No. :	5020-14068
Eluent	: A) CH ₃ CN B) 20 mM KH ₂ PO ₄ (pH 3.0, H ₃ PO ₄) A/B = 25/75 -2.5 min- 25/75 -1.5 min-70/30 -4 min- 70/30 -1 min-25/75 -3 min- 25/75, v/v : 1.0 mL/min
Flow rate	
Column Temp.	: 25 °C
Detection	: UV 240 nm
Injection Vol.	: 5 μL


Tablet Sample

The sample concentration was adjusted to 150 mg/L and 75 mg/L as aspirin (acetylsalicylic acid) and clopidogrel, respectively.



Example of Pretreatment

Method for Tablet Sample



GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences, Inc. Japan

22-1 Nishishinjuku 6-Chome
Shinjuku-ku, Tokyo,
163-1130, Japan
Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsciences.jp
Web: www.glsciences.com

GL Sciences B.V.

De Sleutel 9
5652 AS Eindhoven
The Netherlands
Phone: +31 (0)40 254 95 31
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (ShangHai) Ltd.

Tower B, Room 2003,
Far East International Plaza,
NO,317 Xianxia Road,
Changning District.
Shanghai, China P.C. 200032
Phone: +86 (0)21-6278-2272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn

GL Sciences, Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
Phone: 310-265-4424
Fax: 310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

International Distributors

Visit our Website at:

<https://www.glsciences.com/company/distributor.html>