

# How to Use Preparative HPLC - Part 2


## Scaling up from Analytical HPLC

It is not easy to find out optimal condition for preparative HPLC. Not only large volume of solvent but also substantial amount of precious sample may be required for the evaluation of separation conditions, particularly in preparative HPLC. Consequently, we recommend that the evaluation should be carried out using analytical column (4.6 mm I.D.) in the beginning. Condition for preparative HPLC can be investigated efficiently by using analytical column packed with the same gel as in preparative HPLC column.

In this note, Inertsil ODS-3 was taken as an example, and how to scaling up from analytical column to preparative column is described.

(K. Kanno)

Conditions of preparative HPLC are generally optimized as the following flowchart








<What is important to scale up>

It is important to calculate ratio of cross-sectional area of preparative column to that of analytical column. The ratio can be used as follows;

**1) Increase flow rate in proportion to cross-sectional area of column**

**2) Increase loading amount (injection volume) in proportion to cross-sectional area of column**


					
Inner diameter	4.6 mm	10 mm	14 mm	20 mm	30 mm
Cross-sectional area	16.6 mm <sup>2</sup>	78.5 mm <sup>2</sup>	154 mm <sup>2</sup>	314 mm <sup>2</sup>	707 mm <sup>2</sup>
Ratio of the area <sup>※1</sup>	1	4.7	9.3	18.9	42.5

※1 Ratio of cross-sectional to that of 4.6-mm I.D. column was calculated.

### 1) Increase flow rate in proportion to cross-sectional area of column

The figure shown right indicates relation between linear velocity of mobile phase and height equivalent to theoretical plate (HETP) obtained with three columns packed with 10 μm particles.

Optimum flow rate, at which the lowest HETP is obtained, is 3.0 cm/min (0.5 mm/sec) for all the columns. Therefore, it can be said that flow rate should be changed to maintain optimum linear velocity of 3.0 cm/min in case of scaling up from 10 μm particle packed analytical column to 10 μm particle packed preparative one. It is important that particle size of the two columns is same because optimum flow rate changes also depending on particle size.




### 2) Increase loading amount (injection volume) in proportion to cross-sectional area of column

The figure shown below indicates relation between loading amount and number of theoretical plates (N). Three columns with different inner diameters were used and compared. For example, maximum loading amount for each column at which N above 2000 can be obtained is follows;


- 6 mm I.D. approx. 1 mg**
- 20 mm I.D. approx. 10 mg**
- 50 mm I.D. approx. 70 mg**

Since the maximum loading amount is proportional to cross-sectional area of column, it can be said that similar separation should be achieved with preparative column as with analytical column by increasing injection volume in proportion to cross-sectional area.



<Parameters to be changed for scaling up>

In case of scaling up from a 4.6 mm I.D. analytical column to a 20 mm I.D. preparative column, cross-sectional area of the column is approximately 19 times enlarged. Therefore, scaling up can be achieved by increasing flow rate and loading amount (injection volume) 19 times. **Red letters** represent parameters to be changed for scaling up.

		<b>Scale up</b>		
<b>Column</b>	: Inertsil ODS-3 (10 μm, <b>4.6 mm I.D.</b> × 250 mm)	 Increase flow rate and injection volume 19 times. Keep other parameters unchanged	<b>Column</b>	: Inertsil ODS-3 (10 μm, <b>20 mm I.D.</b> × 250 mm)
<b>Eluent</b>	: A) CH <sub>3</sub> CN B) H <sub>2</sub> O A/B = 40/60, v/v		<b>Eluent</b>	: A) CH <sub>3</sub> CN B) H <sub>2</sub> O A/B = 40/60, v/v
<b>Flow rate</b>	: <b>500 μL/min</b>		<b>Flow rate</b>	: <b>9.5 mL/min</b> (9500 μL/min)
<b>Column Temp.</b>	: 40 °C		<b>Column Temp.</b>	: 40 °C
<b>Detection</b>	: UV 270 nm		<b>Detection</b>	: UV 270 nm
<b>Injection Vol.*</b>	: <b>500 μL</b>	<b>Injection Vol.*</b>	: <b>9.5 mL (9500 μL)</b>	


\* Concentration of the sample solution is same

<An example of scaling up>



Chromatograms before and after scaling up are shown below

**HPLC conditions**

**Column** : Inertsil ODS-3  
(5 μm, 250 mm length)  
**Column I.D.** : Discrised in each chromatogram  
**Eluent** : A) CH<sub>3</sub>CN B) 0.1 % TFA  
A/B = 40/60, v/v  
**Flow rate** : Discrised in each chromatogram  
**Column Temp.** : 40 °C  
**Detection** : UV 270 nm  
**Sample** : Commercially available licorice extract (2.0 mg/mL)



Loading amount was 18.9 times increased.



GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

**GL Sciences, Inc. Japan**

22-1 Nishishinjuku 6-Chome  
 Shinjuku-ku, Tokyo,  
 163-1130, Japan  
 Phone: +81-3-5323-6620  
 Fax: +81-3-5323-6621  
 Email: [world@glsc.co.jp](mailto:world@glsc.co.jp)  
 Web: [www.glsciences.com](http://www.glsciences.com)

**GL Sciences B.V.**

De Sleutel 9  
 5652 AS Eindhoven  
 The Netherlands  
 Phone: +31 (0)40 254 95 31  
 Email: [info@glsciences.eu](mailto:info@glsciences.eu)  
 Web: [www.glsciences.eu](http://www.glsciences.eu)

**GL Sciences (ShangHai) Ltd.**

Tower B, Room 2003,  
 Far East International Plaza,  
 NO,317 Xianxia Road,  
 Changning District.  
 Shanghai, China P.C. 200032  
 Phone: +86 (0)21-6278-2272  
 Email: [contact@glsciences.com.cn](mailto:contact@glsciences.com.cn)  
 Web: [www.glsciences.com.cn](http://www.glsciences.com.cn)

**GL Sciences, Inc. USA**

4733 Torrance Blvd. Suite 255  
 Torrance, CA 90503  
 Phone: 310-265-4424  
 Fax: 310-265-4425  
 Email: [info@glsciencesinc.com](mailto:info@glsciencesinc.com)  
 Web: [www.glsciencesinc.com](http://www.glsciencesinc.com)

**International Distributors**

Visit our Website at:  
<https://www.glsciences.com/company/distributor.html>