LT192 GL Sciences Inc.

Malic acid is one of the organic acids and is widely present in foods and living organisms. Analysis of the malic acid standard reagent confirmed the phenomenon that two peaks appeared depending on the detection method. This time, we analyzed malic acid using a diode array detector, and it was estimated that the impurity contained in malic acid was fumaric acid. (R. Takahashi)

Example. Analysisof standard (1) (Post-column used BTB)



Conditions	
Colum	: InertSphere FA-1
	(9 mm, 300 mm × 7.8 mm I.D.)
Guard Column	: InertSphere FA-1 Guard
	(9 mm, 50 mm × 6.0mm I.D.)
Eluent	: 3 mM HClO ₄ in H ₂ O
Reaction Sol.	: 0.1 mM BTB + 30 mM Na ₂ HPO ₄
Flow Rate	: 0.5 mL/min
Column Temp.	: 35 °C
Detection	: VIS 440 nm
Injection. Vol.	: 50 μL

Example: Analysis of standard (2) (UV absorption)



Conditions

Colum	: InertSphere FA-1
conum	(9 mm, 300 mm × 7.8 mm I.D.)
Guard Column	: InertSphere FA-1 Guard
	(9 mm, 50 mm × 6.0 mm I.D.)
Eluent	: 3mM HClO₄ in H₂O
Flow Rate	: 0.5 mL/min
Column Temp.	: 40 °C
Detection	: UV 210 nm
Injection Vol.	: 10 μL

A peak was observed at UV 280 nm at a position that was barely detected by the BTB post-column method.

>>> Matches the retention time of fumaric acid



UV spectrum of each peak



Malic acid and fumaric acid are similar in structure, but fumaric acid has a conjugated double bond, so the UV absorption wavelength is larger than that of malic acid.

>>> A small amount of fumaric acid is hardly detected by the BTB post-column method, but it is detected as a large peak by the UV absorption method.

HPLC Column

- Analytical Column
 - InertSphere FA-1 9 mm, 300 mm X 7.8 mm I.D. Cat.No. 5020-11003
- Guard Column
 - InertSphere FA-1 Guard 9 mm, 50 mm X 6.0 mm I.D. Cat.No. 5020-10998

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan 22-1 Nishishinjuku 6-chome Shinjuku-ku, Tokyo 163-1130, Japan

Phone: +81-3-5323-6620 Fax: +81-3-5323-6621 Email: <u>world@gls.co.jp</u> Web: www.glsciences.com



International Distributors Visit our Website at www.glsciences.com/distributors

GL Sciences Inc. USA 4733 Torrance Blvd. Suite 255 Torrance, CA 90503 USA

 Phone:
 +1-310-265-4424

 Fax:
 +1-310-265-4425

 Email:
 info@glsciencesinc.com

 Web:
 www.glsciencesinc.com

GL Sciences B.V. Dillenburgstraat 7C 5652AM, Eindhoven The Netherlands

Phone: +31-40-254-9531 Email: info@glsciences.eu Web: www.glsciences.eu

GL Sciences (Shanghai) Limited

Tower B, Room 2003 Far East International Plaza No.317 Xianxia Road, Changning District Shanghai, China 200051

Phone: +86-21-62782272 Email: contact@glsciences.com.cn Web: www.glsciences.com.cn