

サンプリングバッグ外寸 (寸法許容公差 ±5 mm程度)


【フッ素樹脂バッグ】

容量(L)	A辺 (mm)	B辺 (mm)
0.1	85.0	120.0
0.5	125.0	170.0
1	159.0	254.5
2	214.0	254.5
3	254.5	258.0
5	272.5	350.0
10	345.0	460.0
20	475.0	509.0
30	509.0	608.0
50	608.0	700.0
100	708.0	1018.0
200	1204.0	908.0
300	1454.0	1018.0
500	1904.0	1018.0
1000	2500.0	1018.0

【SMART BAG 2F】

2Fは熱溶着による収縮があるため、仕上がり寸法は若干小さくなる。

容量(L)	A辺 (mm)	B辺 (mm)
0.1	84.0	118.0
0.2	84.0	158.0
0.3	105.0	160.0
0.5	130.0	168.0
1	180.0	250.0
2	240.0	250.0
3	300.0	250.0
5	250.0	430.0
10	330.0	500.0
20	510.0	500.0
30	500.0	670.0
50	500.0	920.0
100	750.0	1,000.0
200	1,125.0	1,000.0
300	1,520.0	1,000.0
500	2,330.0	1,000.0



【アルミニウムバッグ】

容量(L)	A辺 (mm)	B辺 (mm)
0.1	84.0	118.0
0.2	84.0	158.0
0.3	105.0	160.0
0.5	130.0	168.0
1	200.0	212.5
2	265.0	212.5
3	212.5	350.0
5	350.0	283.3
10	375.0	425.0
20	425.0	600.0
30	425.0	780.0
50	540.0	850.0
100	875.0	850.0
200	1,370.0	850.0
300	1900.0	850.0
500	3000.0	850.0

【SMART BAG PA】

容量(L)	A辺 (mm)	B辺 (mm)
0.1	84.0	118.0
0.2	84.0	158.0
0.3	105.0	160.0
0.5	130.0	168.0
1	165.0	250.0
2	225.5	250.0
3	270.0	250.0
5	250.0	400.0
10	320.0	500.0
20	495.0	500.0
30	500.0	650.0
50	500.0	900.0
100	735.0	1,000.0
200	1,105.0	1,000.0
300	1,500.0	1,000.0
500	2,300.0	1,000.0

【テドラーバッグ】

容量(L)	A辺 (mm)	B辺 (mm)
0.1	84.0	118.0
0.2	84.0	158.0
0.3	105.0	160.0
0.5	130.0	168.0
1	159.0	254.5
2	214.0	254.5
3	254.5	258.0
5	272.5	350.0
10	345.0	460.0
20	475.0	509.0
30	509.0	608.0
50	608.0	700.0
100	708.0	1,018.0
200	1,204.0	908.0
300	1,454.0	1,018.0
500	1,904.0	1,018.0
1000	4000.0	1018.0
1500	2493.5	1981.2
2000	3325.0	1981.2

【SMART BAG PEs】

容量(L)	A辺 (mm)	B辺 (mm)
1	200.0	200.0
2	220.0	250.0
3	300.0	250.0
5	330.0	330.0
10	350.0	500.0
20	500.0	500.0

【スカイピアバッグ】

容量(L)	A辺 (mm)	B辺 (mm)
1	165.0	250.0
5	250.0	400.0
10	320.0	500.0
50	500.0	900.0