


Volatile organic compounds



System : GC/MS P&T/GC
注) 高速分析用 HTP セットが必要になります。
 Column : AQUATIC-H
 Injection : P&T
 Detection : MS

- Solutes :
1. 1,1-Dichloroethene
 2. Dichloromethane
 3. t-1,2-Dichloroethene
 4. cis-1,2-Dichloroethene
 5. Chloroform
 6. 1,1,1-Trichloroethane
 7. Carbone Tetrachloride
 8. 1,2-Dichloroethane
 9. Benzene
 10. Trichloroethene
 11. 1,2-Dichloropropane
 12. Bromodichloromethane
 13. cis-1,3-Dichloro-1-propene
 14. Toluene
 15. t-1,3-Dichloro-1-propene
 16. 1,1,2-Trichloroethane
 17. Tetrachloroethene
 18. Dibromochloromethane
 19. m-Xylene
 20. p-Xylene
 21. o-Xylene
 22. Bromoform
 23. p-Dichlorobenzene
 - I.S. p-Bromofluorobenzene