

# InertSearch™ for LC

Inertsil® Applications

## Simultaneous analysis of metabolites using LC/MS/MS (5) - compounds related to central metabolism part 2 -

Data No. LL015-0000

*The chromatograms were provided by Yudai Dempo, Takeshi Bamba, and Eiichiro Fukusaki, Department of Biotechnology, Graduate School of Engineering, Osaka University, 2-1 Yamadaoka, Suita, Osaka 565-0871, Japan*

### Conditions

<b>Column</b>	: InertSustain C18 (3 $\mu$ m, 150 x 2.1 mm I.D., Metal-free hardware)
<b>Column Cat. No.</b>	: 5020-00541
<b>Eluent</b>	: A) 10 mM Tributylamine + 15 mM CH <sub>3</sub> COOH in H <sub>2</sub> O B) CH <sub>3</sub> OH
<b>Flow rate</b>	: 0.2 mL/min
<b>Col. Temp.</b>	: 45 °C
<b>Detection</b>	: LC/MS/MS (ESI, Negative, MRM)
<b>Injection Vol.</b>	: 3 $\mu$ L
<b>Sample</b>	: Standard solution (Approx. 5 $\mu$ mol/L each)

Time (min)	A (vol%)	B (vol%)
0	100	0
1	100	0
1.5	85	15
3	85	15
8	50	50
10	0	100
11	0	100
11.5	100	0
17	100	0

Analyte	Precursor ion (m/z)	Product ion (m/z)
Dihydroxyacetone phosphate	169	97
Pantothenic acid	218	88
Nicotinic acid	122	78
Succinic acid	117	73
Carbamoyl phosphate	140	97
Malic acid	133	115
Acetyl phosphate	139	79
$\alpha$ -Ketoglutaric acid	145	101
6-Phosphogluconic acid	275	97
3-Phosphoglyceric acid	185	79
2-Keto-3-deoxy-6-phosphogluconic acid	257	97
Shikimate 3-phosphate	253	97
Citric acid	191	87
Fructose 1,6-bisphosphate	339	97
Phosphoenolpyruvic acid	167	79
Ribulose 1,5-bisphosphate	309	97
4-Hydroxy-3-methyl-but-2-enyl pyrophosphate	261	79
Isocitric acid	191	73
Flavin mononucleotide	455	97
2-Isopropylmalic acid	175	115
5-Phosphoribose 1-pyrophosphate	389	291
Pyrroloquinoline quinone	329	241
10-Camphorsulfonic acid	231	80
Isopentenyl pyrophosphate	245	79
Dimethylallyl pyrophosphate	245	79

# InertSearch™ for LC

Inertsil® Applications

## Simultaneous analysis of metabolites using LC/MS/MS (5) - compounds related to central metabolism part 2 -

